

OPPO STRETCHY GRIP ANKLE WRAP

SKU: OPP240100

PRODUCT DESCRIPTION

FEATURES

- Provides localised support and compression, and avoids unnecessary compression to areas where it is not required.
- Tension is fully adjustable to provide from light support to total compression.
- Cohesive fabric, causing no pain or irritation to the skin whilst in use, comfortable to wear and painless to remove.

HOW TO WEAR

- Place wrap on foot and let the silicone side face to skin.
- Wrap one complete turn around the arch of the foot.
- Angle upward to wrap on complete turn around the lower ankle.
- Continue wrapping in ""figure eight"" pattern.
- Finish and fasten the hook closure.

APPLICATION

- Weak Ankle
- Ankle Fatigue

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au