




OAPL SHORT ANKLE WALKER MOON BOOT

\$102.00



SKU: ORS-

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The OAPL Short Ankle Walker is designed to immobilise the foot and offload forefoot injuries. Featuring Velcro closures, a breathable liner and a short frame with rocker sole to stabilise and reduce ankle movement.

Suitable for metatarsal fractures, plantar fascia tears, bunioneictomies and forefoot soft tissue injuries. Not suitable for 5th metatarsal injuries.

SIZING

Size	Womens (US)	Mens (US)	Footbed Length (cm)
X Small	4 - 5.5	Up to 4	22.5
Small	6 - 8	4.5 - 6.5	26
Medium	8.5 - 10.5	7 - 9	29
Large	11 - 12+	9.5 - 12	30.5
X Large	12.5 - 14	32.5	