



OAPL SHORT ANKLE WALKER MOON BOOT

\$102.00 Original price was: \$102.00.<u>\$92.00</u>Current price is: \$92.00.



SKU: ORS-

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The OAPL Short Ankle Walker is designed to immobilise the foot and offload forefoot injuries. Featuring Velcro closures, a breathable liner and a short frame with rocker sole to stabilise and reduce ankle movement.

Suitable for metatarsal fractures, plantar fascia tears, bunioneictomies and forefoot soft tissue injuries. Not suitable for 5th metatarsal injuries.

SIZING

Size	womens	(US) Mens	(US) Footbed	Lengtn	(cm)
------	--------	-----------	--------------	--------	------

X Small 4 - 5.5 Up to 4 22.5 Small 6 - 8 4.5 - 6.5 26 Medium 8.5 - 10.5 7 - 9 29 Large 11 - 12 + 9.5 - 12 30.5 X Large 12.5 - 14 32.5