




## **BYACRE CARBON FIBRE ULTRALIGHT BACKREST**

**\$99.00**



**SKU: 20105**

## **PRODUCT DESCRIPTION**

## **PURPOSE & BENEFITS**

A backrest that supports your lower back while sitting on the seat of your rollator.

## **FEATURES**

- Great support - for your lower back while sitting on your rollator
- Adjustable height - so it fits no matter the height on your rollator (though limited to a max height 25mm below the rollators tallest position)
- Easy attachment - The backrest is attached into the holes on the handle bar. Put the metal knob on the backrest into one of the holes and then put the Velcro strap on the backrest through the metal bracket and tighten it.

This product is for the [byACRE Ultralight](#) and [byACRE Carbon Overland](#)