

2 AMP SCOOTER BATTERY CHARGER

\$135.00



SKU: TE-CH2



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This is an Australian Approved charger to suit the small Shoprider scooter and Powerchairs (suits other scooters as well). The charger has an LED that glows red whilst charging and green when the Shoprider battery is fully charged.

SPECS

Input 100-240VAC 50-60Hz 2.0A

Output +29.5V Max 60W to suit LITHIUM BATTERY

3 Pin Connector