



LITTMANN STETHOSCOPE CLASSIC III - 3M STETHOSCOPE



SKU: LITTMANN-5620-H

PRODUCT DESCRIPTION

LITTMANN CLASSIC III PURPOSE AND BENEFITS

The Littmann Classic III stethoscope is the most recent model of the stethoscope which helps millions of healthcare professionals achieve their very best.

It brings new layout, materials, and engineering into the Littmann stethoscope series. Using a Littmann Classic III stethoscope illustrates what is important for you as a health professional: grade, quality and performance excellence.

The brand-new Classic III stethoscope provides you with:

- Open bell flexibility - The paediatric side extends to a conventional open bell by substituting the single-piece diaphragm using a non-chill rim
- High Definition sensitivity for doing overall physical examinations
- Tunable diaphragms on either adult or paediatric Surfaces of the chest piece
- Single-piece diaphragm that is simple to attach and easy to clean surface, smooth without any fractures
- Next-generation tubing with enhanced resistance to skin oils and alcohol to get a longer tube lifetime.
- Warranty interval - 5 Years (The 3M Littmann warranty covers the whole device)

The 3M Littmann Classic III stethoscope is created in the USA for constant quality. It comes complete with big and small soft-sealing ear tips, non-cill bell sleeve, and directions.

Stethoscope Characteristics

- Tunable tech is a 3M invention that allows you to change between hearing low and high-frequency sounds simply by adjusting how much you really press the chest piece.
- Soft-sealing ear-tips provide you with a comfortable fit and a superb acoustic seal; they snap tightly on the ear tubes for security.
- Open bell remains clear of debris and dirt from keeping it covered with the little diaphragm
- The flexible headset lets you set the strain for head dimensions and relaxation by simply squeezing or pulling the ear-tubes
- No natural rubber latex or phthalate plasticisers are utilised in the tube or some other element - to help safeguard human health and the environment

Classic III Stethoscope Specifications

Length: 69 cm

Weight: 150g

Applications/occupations: EMT/EMS, Medical Assistant, Nurse, Nursing Student

Acoustic performance: 7

Special adaptors: Yes

Warranty period: 5 Years

CHEST PIECE

Construction: Two-sided, adult and paediatric

Material: Stainless steel

Weight: 82g

Diaphragm type: Tunable, single-piece

Diaphragm diameter: Adult: 4.3cm, Paediatric: 3.3cm

TUBING

Construction: Single lumen (single acoustic channel within tube)

Material: Not made with natural rubber latex or phthalate plasticisers

HEADSET

Material: Wide-diameter aerospace alloy / anodised aluminium

EAR TIPS

Ear tip type: Soft Sealing

Extra Ear tips: Yes

3M™ Littmann® Classic Stethoscopes are constructed to be small but accurate instruments for monitoring and assessing children or adults. They are manufactured with strong yet lightweight materials, allowing you to comfortably pick up subtle heart, lung and other body sounds.

Their soft-sealing ear tips conform to individual ears for an excellent acoustic seal and a comfortable fit. Tunable technology lets you hear high- and low-frequency sounds without needing to reposition your chest piece. The tubing also offers a longer product lifetime by being resistant to skin oils and alcohol and is built with a resilient, latex-free next-generation material that helps retain shape and flexibility.

Buy the Littmann Stethoscope Classic iii with Afterpay or [NDIS method](#) and also, have a look at our other [medical equipment](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au