

## ABLE LIFE UNIVERSAL STANDING HANDLE



**\$38.00**



**SKU:** N/A

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Universal Standing Handle provides a safe way to assist a loved one to stand from a seated position. The 9.5" x 4.5" ribbed plastic handle allows for a comfortable and secure grip for both caregiver and patient. The small standing device is easy to store in a glove box, middle console, or take in a purse to use anywhere. The Universal Stand Assist can support weights of up to 158kg

- **INCREASE SAFETY:** The Universal Standing Handle helps users come to a standing position using maximum leverage and eliminates tears on fragile skin.
- **AVOID INJURIES:** By eliminating the need for under-arm grabbing and pulling with force, the Universal Standing Handle protects against strain on a caretaker's back and prevents harmful pulling on shoulders, arms, and wrists.
- **EASY TO DISINFECT:** Unlike gait belts or lift belts, the handle can be used while keeping a comfortable distance from patients and can easily be sanitized between uses.
- **GREAT FOR CARETAKERS:** The handle is perfect for caretakers, nurses, and therapists to help disabled, handicapped, and senior patients with limited mobility.

### FEATURES

- **INCREASE LEVERAGE:** The Universal Standing Handle provides a safe way to assist a loved one to stand from a seated position.
- **COMFORTABLE GRIP:** The 9.5" x 4.5" ribbed plastic handle allows for a comfortable and secure grip for both caregiver and patient.
- **PORTABLE:** The small standing device is easy to store in a glove box, middle console, or take in a purse to use anywhere.
- **HEAVY DUTY:** The Universal Stand Assist can support weights of up to 350 pounds.
- **THREE COLORS:** Choose between three colour options: Red, Blue, and Ivory.

## **SPECIFICATIONS**

- Weight: 1 lbs
- Weight Capacity: 350 lbs
- Length: 9.5"
- Width: 4.5"