



## ABRI-SOFT CLASSIC UNDERLAY

**\$46.00 - \$70.00** Price range: \$46.00 through \$70.00



**SKU:** ABE4123-H

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

Abri-Soft Classic is a disposable fluff-filled underlay with a waterproof polyethylene backing and sealed edges to ensure leakage control. Made from highly absorbent material, with embossed channels for maximum fluid dispersion and a non-woven layer for optimum dryness and softness. They are suitable for the protection of beds and chairs, blood sampling and changing a baby's nappy when on the move. Products from this range are designed to have minimal impact on the environment and materials used can be recycled.

## FEATURES

- Waterproof backing and sealed edges
- Embossed channels (without flaps) for maximum fluid dispersion
- Non-woven top layer for optimum dryness and softness
- Contains superabsorbent gel (with flaps only)

## SPECIFICATIONS

- 40x60cm 900ml (1pkx60)
- 60x60cm 1300ml (1pkx25)
- 60x90cm 2100ml (1pkx25)