



ABRI-SOFT CLASSIC UNDERLAY

\$41.00 - \$63.00Price range: \$41.00 through \$63.00



SKU: ABE4123-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Abri-Soft Classic is a disposable fluff-filled underlay with a waterproof polyethylene backing and sealed edges to ensure leakage control. Made from highly absorbent material, with embossed channels for maximum fluid dispersion and a non-woven layer for optimum dryness and softness. They are suitable for the protection of beds and chairs, blood sampling and changing a baby's nappy when on the move. Products from this range are designed to have minimal impact on the environment and materials used can be recycled.

FEATURES

- · Waterproof backing and sealed edges
- Embossed channels (without flaps) for maximum fluid dispersion
- Non-woven top layer for optimum dryness and softness
- Contains superabsorbent gel (with flaps only)

SPECIFICATIONS

- 40x60cm 900ml (1pkx60)
- 60x60cm 1300ml (1pkx25)
- 60x90cm 2100ml (1pkx25)