



ACCESS HEALTH DIGITAL INTERMITTENT COMPRESSION UNIT - LEG

\$798.00



SKU: OE1020

PRODUCT DESCRIPTION

Purpose & Benefits

The Access Health Digital Intermittent Compression Unit is designed to support the treatment of lymphoedema, poor circulation, and venous insufficiency by delivering gentle, controlled pneumatic compression to the limbs. This non-invasive therapy enhances blood flow, encourages lymphatic drainage, and helps reduce swelling and discomfort in affected areas.

With a simple digital interface and adjustable pressure settings, the unit is suitable for both clinical and home environments under medical guidance. Its quiet operation and customizable treatment cycles make it an effective and user-friendly solution for managing chronic conditions or aiding in post-operative recovery.

Features

- Digital display with easy-to-use controls
- Adjustable pressure levels and treatment durations
- Promotes blood circulation and lymphatic drainage
- Helps reduce limb swelling and venous pooling
- Suitable for upper and lower limb garments
- Quiet, compact, and portable design

- Ideal for clinics, aged care, and home use
- Manufactured by Access Health

Specifications

- Brand: Access Health
- Power Supply: 220-240V / 50Hz
- Pressure Range: Typically 30-250 mmHg (adjustable)
- Control Panel: Digital with LED indicators
- Timer Function: Adjustable treatment cycles
- Compatibility: Works with standard compression sleeves (sold separately)
- Application: Lymphoedema, venous insufficiency, post-surgical recovery
- Unit Type: Intermittent pneumatic compression (IPC) device