




ACTIMOVE TRI-PANEL KNEE IMMOBILISER

\$87.00Online Only



SKU: 72819

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Actimove Genu Tri-Panel is generally used to reduce pain following trauma or surgery.

The design allows splint adjustments for a broad spectrum of patient leg circumferences and the hook and loop sensitive closure system permits easy application and removal.

FEATURES

- The Actimove Genu Tri-Panel is generally used to reduce pain following trauma or surgery.
- The design allows splint adjustments for a broad spectrum of patient leg circumferences and the hook and loop sensitive closure system permits easy application and removal.
- There are removable and adjustable aluminium stays (two posterior and two on each side panel) which allow a good anatomical fit.
- For patient comfort the Actimove Genu Tri-Panel has a terry cloth interior padding.

SIZING