



ACTION PILOT CUSHION WITH BASIC COVER

\$395.00 - \$491.00 Price range: \$395.00 through \$491.00



SKU: 9000-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ACTION Pilot Cushion is Low-Profile cushion that provides exceptional pressure relief for users with a history of pressure sores. 2.54cm Akton Polymer provides immersion for bony prominences and protects clients with fragile skin integrity.

This Multi-use Cushion is recommended for use on Manual and Power Wheelchairs, or any other seat requiring significant reduced Hi-Pressure areas, vibration and shear.

The Low-Profile allows users to effectively reach floor for more efficient cushions to maximize stability and positioning when required. Custom sizes are available.

The cushion comes with a basic cover, though an incontinent cover is also available.

KEY FEATURES

- The very low profile allows users to effectively reach the floor for more efficient foot-assisted propulsion
- It also conforms to pre-contoured foam or rigid cushions to maximize stability and positioning when required
- Comes with Basic Cover 2.54cm Akton Polymer
- Incontinent and replacement covers are also available

Buy this ACTION Pilot Cushion online at Medimart Australia using [NDIS payment method](#), AfterPay and ZipPay.