



## **A&D ADVANCED PREMIER TALKING PREMIUM BLOOD PRESSURE MONITOR**

**\$135.00**



**SKU: UA-1030T**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The Advanced Premier Talking Blood Pressure Monitor UA-1030T is a fully automatic, voice-assisted monitor designed to support easy and accurate blood pressure tracking at home. Ideal for users with visual impairments or those who prefer audible guidance, it provides spoken instructions and announces results clearly in English.

With features such as irregular heartbeat detection, pressure rating indicators, and a 90-reading memory, this device offers clinically validated performance in a user-friendly format. It's perfect for daily monitoring, with visual and audio feedback to ensure confidence in every reading.

### **Features**

- Voice-guided operation and audio result announcements in English
- Irregular heartbeat detection
- WHO pressure rating indicator for easy interpretation
- Memory storage for 90 readings with date and time
- Cuff fit error and body movement indicators
- One-touch measurement with large digital display

- Clinically validated accuracy
- Includes medium cuff (22–32 cm) and carry case
- Ideal for visually impaired users and home health care

## **Specifications**

- Model: UA-1030T
- Measurement Method: Oscillometric
- Memory: 90 readings with date/time stamp
- Display: Large LCD + voice readout
- Cuff Size: Medium (22–32 cm)
- Power Source: 4 × AA batteries (included) or optional AC adapter
- Audio Output: English language voice guide
- Validation: Clinically tested (BHS/ESH standards)