



## THERA-MED ADJUSTA THERAPEUTIC CUSHION WEDGE



**SKU:** A11503200-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Adjusta Wedge for Body and **Back Support** is a great choice for those seeking comfort management whilst sleeping, in rehabilitation or on bed rest. This product helps elevate your body through its two-part construction and can also be repositioned between 15 degrees and 75 degrees to provide optimum support.

### The Adjusta Wedge can help you:

- Reduce coughing and help drain mucus from the throat to ease symptoms of colds, congestion and post nasal drainage.
- Reduce symptoms of GERD (Gastroesophageal Reflux Disease), Barrett's Oesophagus, Hiatus Hernia, Acid Reflux or Gastric Reflux.
- Assists with breathing difficulties, with those that suffer CHF (Congestive Heart Failure) and other cognitive conditions.
- Can help with snoring and mild sleep apnoea.
- Relieve respiratory discomfort and ease breathing problems due to emphysema, asthma, chest congestion, bronchial, GERD and reflux.
- Alleviate lower back pain through pressure care.

The Adjusta Wedge can also be used with the Knee Relaxer for Back Support, which is designed to release pressure on the legs by lifting the knees into a 'zero-gravity' like position. Both products can be used together or on their own to target areas needing the most support.

## FEATURES

- Allows you to sleep and rest in an elevated position at any angle between 15 degrees and 75 degrees.
- Readily movable, easy to pack away and simple to use.
- Perfect for reading, relaxing, television and recovery from illness or surgery.
- 100% Australian Made.
- Provides good support for head and neck.
- Adjusts up or down as required.
- Can be used in conjunction with the Knee Relaxer for a total zero gravity feel.
- Suitable for use with any contour pillow or neck support.
- Relieves uncomfortable back pain and neck stress when lying flat in bed.
- Supports your body in recommended rest position.
- Adjusta Wedge will keep you elevated and help increase your blood circulation.
- Will reduce acid reflux and help manage GERD while sleeping.

## TWO COVER OPTIONS AVAILABLE

- Quilted Fabric: Attractive, soft, durable and washable polyester fabric.
- SteriPlus Fabric: Moisture and bacteria resistant, features a convenient zip for easy removal, and can be wiped down with a cloth or removed and put in the washing machine.

## SPECIFICATIONS

|  |  |
|--|--|
| <b>Outer Cover SteriPlus</b>                           | Moisture and Bacteria Resistant Polyurethane ply bonded to a knitted polyester backing |
| <b>TGA Listed</b>                                      | Yes  |
| <b>Australian Register of Therapeutic Goods Number</b> | No.174456 Device:Pillow  |
| <b>Length (of back rest portion)</b>                   | 70cm   |
| <b>Width</b>   | 60cm   |
| <b>Size of Semi Circle Support</b>                     | 60cm wide x 70cm high x 18cm deep  |
| <b>Size of Triangular Support</b>                      | 60cm wide x 40cm high x 27cm deep  |
| <b>Inner</b>   | High Density Open Cell Euthane Foam  |
| <b>Outer Cover Quilted Option</b>                      | 100% Cotton quilted cover  |
| <b>Function</b>  | 2 piece construction   |
| <b>Angle Option</b>                                    | 15 degrees to 75 degrees angle   |
| <b>Outer Cover</b>                                     | Tailored Fitted Easy Wash Cotton Cover   |
| <b>Country of Manufacture</b>                          | Australian Made  |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)