


ADJUSTABLE FINGER SPLINT



SKU: PA-10321-H

PRODUCT DESCRIPTION

Purpose & Benefits of Rehacare Adjustable Finger Splint

The Rehacare Adjustable Finger Splint is designed to provide support and stabilization for injured or painful fingers. This splint is particularly beneficial for individuals recovering from fractures, dislocations, sprains, or those dealing with tendon injuries and conditions such as trigger finger or mallet finger. By immobilizing the affected finger, the splint helps to prevent further injury and aids in the proper alignment and healing of the finger, while still allowing some mobility of the hand for daily activities.

Features

- **Adjustable Design:** Features adjustable straps that allow for a custom fit, accommodating different finger sizes and swelling that may occur during the recovery process.
- **Durable Construction:** Made from high-quality, sturdy materials that provide adequate support while remaining comfortable for long-term wear.
- **Padded Interior:** Includes a soft, padded lining that reduces friction against the skin, preventing irritation and enhancing comfort.
- **Ease of Use:** The splint is designed for easy application and removal, which can be done single-handedly, making it convenient for patients with limited dexterity.
- **Versatility:** Suitable for use on any finger, including the thumb, and can be adjusted to fit both the left and right hands.

Specifications

- **Material:** Lightweight yet durable materials, typically featuring a combination of rigid external components and soft internal padding.

- **Size:**

- Size 1 - 2.25 Inches or Less
- Size 2 - 2.25 to 2.5 Inches
- Size 3 - 2.5 to 2.75 Inches
- Size 4 - 2.75 to 3.0 Inches
- Size 5 - 3.0 to 3.25 Inches
- Size 6 - 3.25 to 3.5 Inches

- **Care Instructions:** Easy to clean with a damp cloth; the materials are chosen to ensure durability and hygiene.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au