


ADVANTIFLEX CLASSIC PRESSURE CARE COVER



SKU: ICONSL20MAC-COV

PRODUCT DESCRIPTION

Purpose & Benefits of the Advantiflex Classic Pressure Care Cover

The Advantiflex Classic Pressure Care Cover is a medical-grade mattress cover designed to enhance patient comfort and safety across various healthcare settings. Utilising high-quality materials, this cover provides excellent durability, infection control, and compliance while ensuring optimal patient care.

Features

- **High-Performance Polyurethane Upper Material** – The cover features a single-sided polyurethane upper with 4-way bidirectional stretch, providing shear protection, immersion, and envelopment. It is vapour-permeable, breathable, resistant to patient body fluids, flexible, and gentle on the skin.
- **Heavy-Duty PVC Base Material** – The 100% waterproof base is robust and hard-wearing, protecting the foam core against the bed base. It conforms well to medical bed surfaces, reduces shearing forces, and offers superior resistance to cleaning agents.
- **Infection Control** – All seams are high-frequency sealed and welded, preventing fluid penetration and supporting stringent infection control protocols. The materials are antimicrobial protected, stain-resistant, and can withstand cleaning agents with high-level disinfection compatibility.
- **Compliance** – The cover is fire retardant and meets strict safety standards, making it suitable for healthcare environments requiring high compliance levels.
- **Construction Options** – Available with an integrated zipper with a waterfall protection flap for easy cover removal and fluid ingress protection, or as a fully sealed and enclosed cover without a zipper for high-risk infection control applications.

Specifications

The Advantiflex Classic Pressure Care Cover is a reliable solution for healthcare facilities aiming to provide optimal patient comfort, safety, and infection control.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au