



## AFISCOOTER SE AFIKIM 3 WHEEL SCOOTER



**SKU:** FT00267

## PRODUCT DESCRIPTION

### ABOUT THE AFISCOOTER

The **Afiscooter SE** with its Black & chrome retro design including many features borrowed from the motorcycle industry is an eye-catcher that makes a real statement but its uniqueness doesn't end there, there's plenty of substance too. In addition to the style it brings, the strong steel chassis and strong power unit creates an extremely rugged scooter with a 20 year track record of unbeatable reliability.

The **Afiscooter SE** has the comfort and stability of a big scooter, simple controls and smooth steering that together with an extra small turning radius allows the driver to negotiate almost any obstacle, be it in the high streets or on a countryside traverse. Go to the beach or the country club, take a tour around your village, go to the shopping center, with the Afiscooter Se you can do it all!

### FEATURES

- Ergonomic Tiller with a Reliable & simple electronic command system for easy & smooth operation
- A Luxurious orthopedic rotatable seat with supportive armrest for maximum comfort and ease of access
- Solid large rear basket
- High Ground Clearance together with Big size Pneumatic tires give excellent outdoor capabilities
- 3 wheel model for maximum maneuverability

### SPECIFICATIONS

- Overall Length - 158cm (62")
- Overall Width - Standard Wheels - 73cm (29")
- Overall Height - No Canopy - 124cm (49")
- Overall Height - with Canopy - 202cm (79")
- Ground Clearance - 125cm (5")

- Weight - No Batteries - 99kg (218lb)
- Weight - (+) Batteries - +146kg (322lb)
- Weight Capacity Including Driver - 200kg (550lb)
- Travel Range - 45km (28m)
- Maximum Speed - 15Km/h (9.3mph)\*
- Max Climbing Angle - 20%
- Turning Radius - 115cm (45")
- Seating Width - 46cm (18")
- Tire Type - Pneumatic (1.7bar)
- Tire Size - Front - 2.75x10"
- Tire Size - Rear - 2.75x12"
- Battery Type x2 - 12V 73Ah
- Rear Wheel Drive Power Unit - 24V DC motor | 1,300W
- Charger Type - 8 AMP

Buy This Afiscooter SE Afikim 3 Wheel Scooter at MediMart using Afterpay or [NDIS method](#) today!  
Expert [mobility equipment suppliers](#) in Australia.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)