



AFISCOOTER SE AFIKIM 3 WHEEL SCOOTER

\$7,476.00 Original price was:

\$7,476.00.<u>\$6,728.40</u>Current price is: \$6,728.40.



SKU: FT00267

PRODUCT DESCRIPTION

ABOUT THE AFISCOOTER

The **Afiscooter SE** with its Black & chrome retro design including many features borrowed from the motorcycle industry is an eye-catcher that makes a real statement but its uniqueness doesn't end there, there's plenty of substance too. In addition to the style it brings, the strong steel chassis and strong power unit creates an extremely rugged scooter with a 20 year track record of un beatable reliability.

The **Afiscooter SE** has the comfort and stability of a big scooter, simple controls and smooth steering that together with an extra small turning radius allows the driver to negotiate almost any obstacle, be it in the high streets or on a countryside traverse. Go to the beach or the country club, take a tour around your village, go to the shopping center, with the Afiscooter Se you can do it all!

FEATURES

- Ergonomic Tiller with a Reliable & simple electronic command system for easy & smooth operation
- A Luxurious orthopedic rotatable seat with supportive armrest for maximum comfort and ease of access
- Solid large rear basket
- High Ground Clearance together with Big size Pneumatic tires give excellent outdoor capabilities
- 3 wheel model for maximum maneuverability

SPECIFICATIONS

- Overall Length 158cm (62")
- Overall Width Standard Wheels 73cm (29")
- Overall Height No Canopy 124cm (49")



- Overall Height with Canopy 202cm (79")
- Ground Clearance 125cm (5")
- Weight No Batteries 99kg (218lb)
- Weight (+) Batteries +146kg (322lb)
- Weight Capacity Including Driver 200kg (550lb)
- Travel Range 45km (28m)
- Maximum Speed 15Km/h (9.3mph)*
- Max Climbing Angle 20%
- Turning Radius 115cm (45")
- Seating Width 46cm (18")
- Tire Type Pneumatic (1.7bar)
- Tire Size Front 2.75x10"
- Tire Size Rear 2.75x12"
- Battery Type x2 12V 73Ah
- Rear Wheel Drive Power Unit 24V DC motor | 1,300W
- Charger Type 8 AMP

Buy This Afiscooter SE Afikim 3 Wheel Scooter at MediMart using Afterpay or <u>NDIS method</u> today! Expert <u>mobility equipment suppliers</u> in Australia.