

AFIKIM FEATHERWEIGHT WHEELCHAIR

\$1,329.00



SKU: FW1



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Afikim Featherweight Wheelchair is one of the lightest wheelchairs on the market. It is made with high quality materials that makes it easy to propel and push.

This featherweight wheelchair is a perfect travelling companion and makes it a breeze to get in and out of the car. It can be folded down into a compact package by folding the backrest down and then pulling the middle of the seat up.

Optional quick-release wheels allow a caregiver to remove the larger rear wheels with the push of the button. Once you arrive at your destination, unfolding and attaching the wheels takes less than a minute.

Standard features include wheel locks at the front of the frame which can be easily reached by the user. Wheel locks prevent the wheelchair from rolling forward or backward.

Additionally, there is an integrated braking system operated by bicycle style locking brake levers located at the push handles for the carer.

These handbrakes provide the caregiver control, security, and convenience when descending or pausing on a slope.

FEATURES

- At 10.3kg this is one of the lightest wheelchairs in the world
- Folds tightly for easy travel

- Carry handles
- One-button, quick release rear
- Wheel locks
- Grey seat cushions cover options available
- Optional carry bag
- Full length fix armrests
- Wheel Material - Featherweight Solid, Polyurethane
- Push to lock Brakes
- Removable and Swing Away Elevating Leg rest pops in at the top and locks with the black plastic knob below.
- Standard features include wheel locks at the front of the frame. which can be easily reached by the user. Wheel locks prevent the wheelchair from rolling forward or backward. Additionally, there is an integrated when descending or pausing on a slope.

SPECIFICATIONS

Specification	Measurement
Product Weight (With Rear Wheels)	10.3kg
Product Weight (Without Rear Wheels)	6.7kg
Maximum Weight Capacity	108kg
Unfolded Size (L x W x H)	76.2 x 70 x 89 cm (30 x 27.6 x 35 in)
Folded Size w/ cushion removed (L x W x H)	58.4 x 35.5 x 73.7 cm (23 x 14 x 29 in)
Frame Material	Aluminium
Travel Distance	Self Propel
Seat Size (Width x Depth)	45.7 x 40.6 cm (18 x 16 in)
Seat to Floor Height w/ cushion	44.5 cm (17.5 in)
Seat to Floor Height w/o cushion	43.2 cm (17 in)
Armrest Type	Fixed
Armrest Height	67.3 cm (26.5 in)
Back height	40.6 cm (16 in)
Turn Radius	80 cm (31.5 in)
Front wheel size	15.2 cm (6 in)
Rear wheel size	55.9 cm (22 in)
Wheel Material	Feather Chair Solid, Polyurethane
Upholstery	Nylon
Footrest	Swing away