




AFIKIM ULTRA LIGHT 18" SELF PROPELLED WHEELCHAIR



SKU: 111431

PRODUCT DESCRIPTION

Purpose & Benefits

The Afikim Ultra Light 18" self-propelled wheelchair is designed to deliver reliable mobility in a lightweight, easy-to-use format. Ideal for users who value independence, this wheelchair features a durable yet lightweight aluminium frame, allowing for easy manoeuvrability both indoors and outdoors.

Its foldable design and quick-release wheels make transport and storage simple, while the padded seat and armrests ensure lasting comfort during use. Whether used at home, in care facilities, or while travelling, it provides a convenient, user-friendly solution for everyday mobility.

Maintain freedom of movement with an ultra-light self-propelled design

Easily folds for car transport or compact storage

Quick-release wheels simplify handling and reduce bulk

Comfortable padded armrests and seat for extended use

Features

- Quick-release rear wheels
- Compact folding design for travel and storage
- Built-in carry handles

- Push-to-lock braking system
- Forward-mounted wheel locks for enhanced safety
- Swing up full-length padded armrests
- Removable and swing-away leg rest
- Optional carry bag available
- Designed for safe use on flat or sloped surfaces
- Lightweight and easy to manage for both users and carers

Specifications

- Seat Width: 18"
- Total Weight: 10.3kg
- Weight Without Rear Wheels: 6.7kg
- Wheel Material: Solid polyurethane featherweight
- Frame: Lightweight aluminium
- Brakes: Push-to-lock and slope-assist wheel locks
- Leg Rest: Removable, swing-away, with top pop-in and lock system
- Maximum User Weight: 115kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au