

AFISCOOTER C 4 WHEEL SCOOTER



SKU: N/A



PRODUCT DESCRIPTION

ABOUT THE AFISCOOTER C

The **Afiscooter C** with its unique design made to maximize the protection and safety of the driver, its excellent multifunctionality and ease of operation is the perfect mobility scooter for those who wish to maintain an active, dynamic and versatile lifestyle.

The combination of Mid-size dimensions, a powerful motor, a luxurious captain seat, well-tuned shock absorbers and an ergonomic modern steering tiller creates a scooter that is brilliant for both Outdoor and Indoors purposes. The three-wheel model has an extremely small turning radius and give supreme maneuverability

Whichever you choose you will enjoy an exceptionally comfortable ride, simple operation and a premium look and feel.

KEY FEATURES

- **User Friendly**, with an ergonomic tiller, left and right-hand controls and an easy to read display
- **Great storage space**, with Built in cup holders, a lockable front compartment and an option for a Rear lock box
- **Multi Purpose**, with Mid-size dimensions, Pneumatic tires, a powerful motor and High Ground Clearance that are brilliant for both Outdoor and Indoors purposes.
- **Safe and Comfortable**, with well-tuned shock absorbers, a Luxurious orthopedic rotatable seat and a Powerful wide field of view LED headlight
- **Excellent Travel Range**, with capacity to carry batteries up to 85 Ah

[Link to Brochure -](#)

https://afiscooters.com/wp-content/uploads/2020/02/AfiScooter-C_Product-Page-2020_4Web.pdf

[Link to User Manuel -](#)

https://afiscooters.com/wp-content/uploads/2020/02/afiscooter-c_-user-manual.pdf

Buy This Afiscooter C 4 Wheel Scooter online at MediMart using Afterpay or [NDIS method](#) today!
Expert [mobility equipment suppliers](#) in Australia.

SPECIFICATIONS

- Available Colours: Silver / Gray / Blue / Red
- Dimensions: 137 x 67 x 80 cm (L x W x H)
- Height Overall Max.: 128 cm
- Height Overall Max. with Canopy: 163 cm
- Weight (total): 141 kg
- Weight without Batteries: 95 kg
- Weight capacity: 150 kg
- Turning radius: 170 cm
- Maximum speed: 10 km/h
- Maximum range *: 40 – 45 km
- Maximum slope *: 9 ° / 16%
- Ground Clearance: 10 cm
- Seating: 46 cm Luxurious Captain style, fully adjustable
- Motor / Gearbox: 24-volt 950 Watt 4 pole DC, sealed transaxle
- Control operation: Wigwag paddle, variable speed, two speed, Delta handle bar
- Dual Braking: Electronic, regenerative and electromechanical
- Electronics: New LCD digital screen
- Tyres: Type: Pneumatic
 - Front : 4.10/3.50-6
 - Rear : 4.10/3.50-6
- Suspension: Full suspension, four points
- Batteries: 2 x 12 V x 55 Ah GEL/AGM Batteries.
- Battery charger: Off board 8 Amp

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au