



AFISCOOTER C4 SCOOTER

\$4,840.00 - \$4,841.88 Price range: \$4,840.00 through \$4,841.88



SKU: FTC4536-B-H

PRODUCT DESCRIPTION

Purpose & Benefits of Afiscooter C4 Scooter

The **Afiscooter C4 Scooter** is a premium 4-wheel mobility scooter designed to deliver a smooth, stable, and highly comfortable ride for users seeking performance and reliability. Built for both urban and rural environments, the C4 is ideal for those who want to maintain an active, independent lifestyle with the support of a robust and stylish scooter.

With full suspension, a luxurious orthopaedic seat, and advanced safety features, the C4 offers an exceptional riding experience on a variety of terrains. Its ergonomic controls and powerful motor ensure confidence and ease of use, whether you're running errands, visiting friends, or simply enjoying the outdoors.

Key Features:

- 4-wheel design for maximum stability on various surfaces
- Full suspension system for a smooth, comfortable ride
- Orthopaedic high-back seat with swivel and sliding adjustments
- Powerful motor and long-range battery for extended travel
- Delta tiller with ergonomic, easy-grip controls

- Full lighting package including indicators and hazard lights for added safety
- Large front basket and optional rear storage box
- Suitable for everyday use in urban and rural settings

Specifications:

- **Maximum Speed:** up to 10 km/hr
- **Front Wheel:** 12 inch
- **Rear Wheel:** 12 inch
- **Dim.:** 1320(L)*670(W)*1260(H) mm
- **Weight Capacity:** 150 kg
- **Per Charge Range:** 40 km
- **Chassis Height Above Ground:** 100 mm
- **Net Weight:** 122 kg
- **Battery Weight (One Unit):** 12V/75AH 15-20kg