


## AFISCOOTER C4 SCOOTER



**SKU:** FTC4536-B-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of Afisooter C4 Scooter

The **Afisooter C4 Scooter** is a premium 4-wheel mobility scooter designed to deliver a smooth, stable, and highly comfortable ride for users seeking performance and reliability. Built for both urban and rural environments, the C4 is ideal for those who want to maintain an active, independent lifestyle with the support of a robust and stylish scooter.

With full suspension, a luxurious orthopaedic seat, and advanced safety features, the C4 offers an exceptional riding experience on a variety of terrains. Its ergonomic controls and powerful motor ensure confidence and ease of use, whether you're running errands, visiting friends, or simply enjoying the outdoors.

### Key Features:

- 4-wheel design for maximum stability on various surfaces
- Full suspension system for a smooth, comfortable ride
- Orthopaedic high-back seat with swivel and sliding adjustments
- Powerful motor and long-range battery for extended travel
- Delta tiller with ergonomic, easy-grip controls

- Full lighting package including indicators and hazard lights for added safety
- Large front basket and optional rear storage box
- Suitable for everyday use in urban and rural settings

## Specifications:

- **Maximum Speed:** up to 10 km/hr
- **Front Wheel:** 12 inch
- **Rear Wheel:** 12 inch
- **Dim.:** 1320(L)\*670(W)\*1260(H) mm
- **Weight Capacity:** 150 kg
- **Per Charge Range:** 40 km
- **Chassis Height Above Ground:** 100 mm
- **Net Weight:** 122 kg
- **Battery Weight (One Unit):** 12V/75AH 15-20kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)