




## AFISCOOTER S4 SCOOTER



**SKU:** FTS4530-B-H

## PRODUCT DESCRIPTION

The **Afiscooter S4 Scooter** is a top-tier, all-terrain mobility scooter built for users who demand the highest level of comfort, performance, and safety. Designed for long-distance travel and outdoor exploration, the S4 offers a smooth, stable ride across various terrains — from city footpaths to country tracks — making it ideal for those who live an active and independent lifestyle.

Equipped with full suspension, an extra-wide orthopaedic seat, and an advanced lighting system, the S4 delivers a luxury driving experience without compromising safety. Its powerful motor and extended battery range ensure reliability for full-day outings, while the robust four-wheel design enhances stability and user confidence in all conditions.

## Key Features:

- Heavy-duty 4-wheel scooter ideal for long-range outdoor travel
- Full suspension system for superior comfort over rough terrain
- Extra-wide, fully adjustable orthopaedic seat with swivel and recline
- Powerful motor with an extended battery range
- Ergonomic delta tiller with easy-to-use controls
- LED lighting package, indicators, and rearview mirrors for safety

- High ground clearance and large pneumatic tyres for off-road capability
- Front and rear storage options for added convenience
- Weight capacity up to 200 kg

## Specifications:

- **Maximum Speed:** 10 km
- **Front Wheel:** 30.48 cm
- **Rear Wheel:** 38.1 cm
- **Dimensions:** 154(L)\*68(W)\*90(L) cm
- **Weight Capacity:** 200 kg
- **Per Charge Range:** 50km
- **Battery:** 12V/75AH (Lead-acid battery)\*2 (Standard)
- **Turning Radius (mm):** 228 cm
- **Chassis Height Above Ground:** 12 cm
- **Seat Height Above Ground:** 66 cm
- **Net Weight:** 155 kg
- **12V/75AH Battery Weight (One Unit):** 15-20kg
- **Seat Weight:** 20kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)