




### AFISCOOTER S4 SCOOTER

**\$7,793.91 - \$8,252.77** Price range: \$7,793.91 through \$8,252.77



**SKU:** FTS4530-B-H

## PRODUCT DESCRIPTION

The **Afiscooter S4 Scooter** is a top-tier, all-terrain mobility scooter built for users who demand the highest level of comfort, performance, and safety. Designed for long-distance travel and outdoor exploration, the S4 offers a smooth, stable ride across various terrains — from city footpaths to country tracks — making it ideal for those who live an active and independent lifestyle.

Equipped with full suspension, an extra-wide orthopaedic seat, and an advanced lighting system, the S4 delivers a luxury driving experience without compromising safety. Its powerful motor and extended battery range ensure reliability for full-day outings, while the robust four-wheel design enhances stability and user confidence in all conditions.

### Key Features:

- Heavy-duty 4-wheel scooter ideal for long-range outdoor travel
- Full suspension system for superior comfort over rough terrain
- Extra-wide, fully adjustable orthopaedic seat with swivel and recline
- Powerful motor with an extended battery range
- Ergonomic delta tiller with easy-to-use controls
- LED lighting package, indicators, and rearview mirrors for safety

- High ground clearance and large pneumatic tyres for off-road capability
- Front and rear storage options for added convenience
- Weight capacity up to 200 kg

## Specifications:

- **Maximum Speed:** 10 km
- **Front Wheel:** 30.48 cm
- **Rear Wheel:** 38.1 cm
- **Dimensions:** 154(L)\*68(W)\*90(L) cm
- **Weight Capacity:** 200 kg
- **Per Charge Range:** 50km
- **Battery:** 12V/75AH (Lead-acid battery)\*2 (Standard)
- **Turning Radius (mm):** 228 cm
- **Chassis Height Above Ground:** 12 cm
- **Seat Height Above Ground:** 66 cm
- **Net Weight:** 155 kg
- **12V/75AH Battery Weight (One Unit):** 15-20kg
- **Seat Weight:** 20kg