




## **AIDACARE ASPIRE A150F ALUMINIUM FOLDING LIFTER**

**\$4,882.00 Original price was:**

**~~\$4,882.00~~ ~~\$4,005.00~~ Current price is: \$4,005.00. Sale Available Online Only!**



**SKU: LSS390600**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Aidacare Aspire A150F Aluminium Folding Lifter**

The **Aidacare Aspire A150F Aluminium Folding Lifter** is a lightweight, compact patient lifter designed to assist with safe and efficient transfers in home care, aged care, and clinical environments. Ideal for lifting and repositioning individuals with limited mobility, the A150F enhances both user safety and carer convenience.

Constructed from durable aluminium, this lifter features a folding frame for easy transport and storage, making it ideal for environments where space is limited, or mobility equipment needs to be relocated frequently. The electric lift mechanism, adjustable leg base, and ergonomic handset allow for smooth, controlled transfers with minimal strain on the caregiver.

### **Key Features:**

- Folding aluminium frame for compact storage and transport
- Electric lift mechanism with ergonomic handset
- Adjustable leg base for easy access around furniture
- Safe working load of 150 kg
- Quick-release leg and mast locking system for fast assembly/disassembly

- Durable castors for smooth manoeuvrability
- Compatible with a wide range of slings
- Ideal for home care, residential aged care, and rehabilitation settings

## Specifications:

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Safe Working Load	Material
LSS390600	Aidacare							
	Aspire							
	A150F	1220 -	660 -	1873	428 -	580 -	150 kg	Aluminium
	Folding Patient Lifter	1280 mm	1040 mm	mm	1640 mm	950 mm		