




AIR COMFORT COMPACT MOBILE LIFT CHAIR DUAL MOTOR



SKU: AC5912-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Air Comfort Compact Mobile Lift Chair – is a harmonious blend of luxury and functionality, meticulously crafted for those who seek exceptional comfort and impeccable design. With its unique features and state-of-the-art technology, this chair is not just a piece of furniture but an experience that redefines relaxation. Whether you're indulging in a leisurely afternoon or require specific support and mobility, this chair caters to all your needs, combining style, convenience, and unmatched comfort. Dive into its many features and elevate your lounging experience to unparalleled heights.

FEATURES

- **Pressure Relief:** At the heart of this chair is the exclusive Air Comfort Seating system. It's engineered to provide long-term pressure relief, ensuring maximum comfort for extended periods.
- **Modern & Easy Maintenance:** Outfitted with the VersoTela Theme material, the chair not only sports a contemporary look that complements any room décor but is also a breeze to clean.
- **Customizable Comfort:** Catering to your specific needs, it has a removable seat and back cushions. Tailor the chair's cushioning to your liking and enjoy customized coziness.
- **Power-Assisted Movements:** With its assisted electric sit, lift, and recline function, transitioning between positions has never been smoother. The twin-motor system ensures efficient and seamless operation, making every movement effortless.
- **Mobility & Positioning:** Fitted with 100mm front and rear castors, you can easily position the chair wherever you want. Whether you're looking to shift it to a cozy corner or bring it centre stage, maneuvering is a cinch.
- **Wall Saver Mechanism:** No more worries about scraping or damaging your walls. The chair's



Wall Saver mechanism ensures it reclines without coming into contact with your wall, preserving both the chair and your interiors.

- **Always Powered:** Unexpected power outages? Not a problem. The backup battery capability ensures the chair's functions aren't interrupted, giving you peace of mind.
- **Durability & Easy Upkeep:** This chair is not only easy to clean but also simple to repair and maintain. Its robust build ensures longevity, and its design means any maintenance is hassle-free.

SPECIFICATIONS

	Dimension	Value
Overall Length		820 mm
Overall Width		750 mm
Overall Height		1060 mm
Seat Width		500 mm
Seat Depth		410 mm
Floor To Seat Height		480 mm
Backrest Height		650 mm
Max User Weight		140 kg
Brand		Air Comfort

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au