



AIR COMFORT COMPACT PRESSURE RELIEF RECLINER LIFT CHAIR

\$3,606.00 - \$3,854.00 Price range: \$3,606.00 through \$3,854.00



SKU: AC590-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience luxurious comfort with this lift and recline chair, which is available in three sizes.

The Air Comfort seating system is designed to aid in the treatment and prevention of pressure problems related to long term seating. Utilising 4 pressure-less, air filled, high frequency welded PVC bags, Air Comfort products work to distribute the occupants weight, thereby ensuring there are no main pressure points restricting circulation.

Assists in the aid of pressure related ailments and is perfect for people with back pain, hip problems, osteo related problems, circulation ailments, degenerative diseases (MS).

FEATURES

- Removable seat and back cushions
- Heavy-duty front castors and rear locking castors
- Cover is fusion Carrflex 4-way stretch Vinyl
- Wall saver mechanism
- Wooden hand knuckles
- Back-up battery capability
- Weight capacity of 140kg
- Twin motor also available for full bed recline, and independent control of back and leg supports.
- Total surface air pressure system
- Electric lift and recline mechanism

FEATURES OF FABRIC

- Antifungal
- Fire resistant
- Anti-bacterial
- Wear-resistant

SPECIFICATIONS

- Small
 - Seat Height - 48cm
 - Seat Width - 50cm
 - Seat Depth - 44cm
 - Seat back Height - 68cm
- Medium
 - Seat Height - 48cm
 - Seat Width - 50cm
 - Seat Depth - 51cm
 - Seat back Height - 74cm
- Large
 - Seat Height - 52cm
 - Seat Width - 60cm
 - Seat Depth - 51cm
 - Seat back Height - 80cm

<https://youtu.be/prNSdedtqxA>