


AIR COMFORT DELUXE BED V2 - ELECTRIC



SKU: AC59112

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Air Comfort Deluxe Bed V2 - Electric features a Patented Airgonomic Seating System that uses four air cells to provide immersive and cushioned comfort and is designed for patients who are between 1620mm to 1880mm in height.

The bed also includes an electric adjustable back recline, seat tilt, and raising legrest, as well as drop-away and adjustable armrests and back wings. The VersoTela Fusion material used in the design is specifically intended for pressure care and sensitive skin.

Manoeuvring the bed is made easy with its large castors, and it includes a one-directional lock castor to help with straight-line tracking. Customisation and colour options are also available.

This is also available in a standard [Air Comfort Deluxe Bed V2](#) model.

FEATURES

- Patented Airgonomic Seating System that uses four air cells for immersive and cushioned comfort
- Electrically adjustable back recline, seat tilt, and raising legrest for personalised positioning
- Drop-away and adjustable armrests and back wings for extra support and comfort
- Features VersoTela Fusion material designed specifically for pressure care and sensitive skin
- Large castors for easy maneuverability and mobility
- One directional lock castor to assist with straight line tracking and stability
- Customisation and colour options available for personalised style and design
- Recommended for patient height up to 1620mm to 1880mm for optimal comfort and support

- Maximum user weight: 180 kg

COLOURS

If you prefer an alternative colour, please enter your preference in the order notes section of the checkout.

ACCESSORIES

The following accessories are compatible with the Air Comfort Deluxe Bed V2 – Maxi. These optional add-ons enhance comfort, safety, and functionality:

- [Footrest for Deluxe Bed V2](#) – Provides lower limb support and pressure care.
- [Meal Tray for Deluxe Bed V2](#) – Offers a sturdy surface for meals, reading, and more.
- [Incontinence Cover for Deluxe Bed V2](#) – Waterproof and breathable protection for the mattress.
- [Seat Belt for Deluxe Bed V2](#) – Improves user safety.

SPECIFICATIONS

Recommended for user height	1620 to 1880mm
Overall Length	1130 mm
Overall Width	680 mm
Overall Height	1260 mm
Seat Width	470 mm
Seat Depth	550 mm
Floor To Seat Height	570 mm
Product Weight	45 kg
Brand	Air Comfort

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au