




AIR COMFORT DELUXE BED V2 - MAXI ELECTRIC



SKU: AC59115

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Air Comfort Deluxe Bed V2 - Maxi is an advanced therapeutic bed designed to provide exceptional comfort and support for individuals requiring extended bed rest or those with specific healthcare needs. This bed combines cutting-edge technology with user-centric design to offer an optimal sleeping and resting environment, particularly for bariatric patients or those at risk of pressure injuries.

FEATURES

- **Advanced Pressure Redistribution:** The bed is equipped with a sophisticated air system that adjusts pressure distribution dynamically, catering to the individual's body contours and movements to minimize pressure points and the risk of pressure ulcers.
- **Customizable Support:** Features multiple air cells that can be individually adjusted to achieve the perfect balance of support and comfort, ensuring that patients with varying needs are adequately accommodated.
- **Integrated Pump System:** A built-in, low-noise pump that continuously regulates airflow and pressure within the mattress and maintains a consistent and comfortable sleeping surface.
- **Easy-to-Use Control Panel:** The user-friendly control panel allows caregivers or users to easily adjust settings, including firmness levels and cycle times, for customized care.
- **Waterproof, Breathable Cover:** Includes a high-quality, stretchable cover that is waterproof yet breathable, promoting air circulation to keep the skin cool and dry, while also being easy to clean and disinfect.
- **Bariatric Capacity:** Specifically designed to support higher weight capacities, making it suitable for bariatric patients while ensuring stability and safety.

ACCESSORIES

The following accessories are compatible with the Air Comfort Deluxe Bed V2 – Maxi. These optional add-ons enhance comfort, safety, and functionality:

- [Footrest for Deluxe Bed V2](#) – Provides lower limb support and pressure care.
- [Meal Tray for Deluxe Bed V2](#) – Offers a sturdy surface for meals, reading, and more.
- [Incontinence Cover for Deluxe Bed V2](#) – Waterproof and breathable protection for the mattress.
- [Seat Belt for Deluxe Bed V2](#) – Improves user safety.

SPECIFICATIONS

Overall Length	1360 mm
Overall Width	950 mm
Overall Depth	950 mm
Overall Height	1385 mm
Seat Width	720 mm
Seat Depth	550 mm
Floor To Seat Height	640 mm
Product Weight	100 kg
Max User Weight	300 kg
Brand	Air Comfort

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au