



AIRCOMFORT DELUXE BED STANDARD - ELECTRIC - INCLUDING POWER DRIVE



SKU: AC59112PD

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aircomfort Deluxe Bed Standard - Electric - Including Power Drive** is designed to provide advanced pressure care and comfort for patients. Using its patented Airgonomic Seating System, it immerses and cushions the user, while powered adjustments and specialised materials support safe positioning, easier transfers and improved skin protection.

Features

- Patented Airgonomic Seating System with four air cells for immersion and cushioning
- Electric adjustments for back recline, seat tilt and raising legrest
- Drop-away and adjustable armrests and back wings for easier transfers
- VersoTela Fusion material designed for pressure care and sensitive skin
- Large castors for easier manoeuvrability
- One directional lock castor assists with straight line tracking
- Customisation and colour options available

Specifications

- **Drive Direction:** Forwards & backwards
- **Max Assist Speed:** Up to 6 km/h
- **Max Advised Ramp Angle:** 6° Wheel Lift Speed: 4 seconds up/down
- **Product Weight (excl. charger):** 18 kg
- **Force of Wheel on Floor:** 40 kg
- **Battery Voltage / Capacity:** 24V DC / 9Ah NiMH
- **Battery Lifetime:** 500 cycles or ~2 years
- **Average Running Distance:** ~9 km on full charge (varies by load, slope & use)
- **Safety Features:** Emergency stop button, overcharge management, IPX6 charger

Key Features

- Powered drive assistance for effortless patient transport
- Enhances workplace safety by reducing manual handling
- Intuitive force-sensing handles for proportional propulsion
- Retractable drive mechanism for manual push option
- Simple user interface with battery indicator, service alert & emergency stop
- IPX6 water resistant battery charger
- Designed exclusively for trained healthcare staff in clinical facilities