



## **AIR COMFORT DELUXE BED V2**

**\$4,084.00**



**SKU: AC592**

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Air Comfort Deluxe Chair V2 is a seating solution designed to provide pressure relief and redistribution. It is a popular choice in the market as it functions as a hybrid day chair and bed, offering maximum comfort for sensitive patients and ease of use for caregivers.

The chair features adjustable back recline, seat tilt, and raising leg rest, as well as drop away/adjustable armrests and back wings to enhance user comfort. Its Carrflex Extreme multi-stretch covering has been treated for sanitation and durability.

The chair is also designed for mobility, with large castors and a one-directional lock castor to facilitate straight-line tracking. The patented Airgonomic seating system of the chair comprises four air cells that provide a cushioned and supportive position for the patient.

Available in three sizes, please choose from the dropdown menu above when placing an order.

This is also available in an [Air Comfort Deluxe Bed V2 - Electric](#) model.

### **FEATURES**

- Recommended patient height: Up to 1620mm to 1880mm.
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the user.
- Adjustable back recline, seat tilt and raising legrest.
- Drop away / adjustable arm rests and back wings.
- Features VersoTela Fusion material specifically designed for pressure care and sensitive skin.
- Large castors make the chair easy to manoeuvre.

- One directional lock castor to assist with straight line tracking.
- Customisation and colour options available.

## COLOUR

If you would prefer an alternative colour, please specify in the order notes section at checkout.

## ACCESSORIES

The following accessories are compatible with the Air Comfort Deluxe Bed V2 – Maxi. These optional add-ons enhance comfort, safety, and functionality:

- [Footrest for Deluxe Bed V2](#) – Provides lower limb support and pressure care.
- [Meal Tray for Deluxe Bed V2](#) – Offers a sturdy surface for meals, reading, and more.
- [Incontinence Cover for Deluxe Bed V2](#) – Waterproof and breathable protection for the mattress.
- [Seat Belt for Deluxe Bed V2](#) – Improves user safety.

## SPECIFICATIONS

Product Name	SMALL	MEDIUM	STANDARD
Overall Length (mm)	1000	1100	1130
Overall Width (mm)	670	670	680
Overall Height (mm)	1130	1130	1260
Seat Width (mm)	430	430	470
Seat Depth (mm)	400	470	550
Floor To Seat Height (mm)	550	550	570
Product Weight (kg)	40	40	45
Recommended patient height (mm)	Up to 1420	1420 to 1620	1620 to 1880
Safe Working Load (kg)	180	180	180