



AIRGO COMFORT-PLUS FOLDING CANE - BLUE

~~\$50.00~~ Original price was: \$50.00. \$45.00 Current price is: \$45.00.



SKU: 730-490

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Airgo Comfort-Plus Folding Cane** is engineered for maximum comfort, security, and convenience. Designed for individuals seeking a lightweight and portable mobility aid, it features an ergonomic cushion-top handle that reduces shock and strain on the hand and wrist. The folding design allows for compact storage and portability, making it ideal for travel or everyday use.

FEATURES

- Ergonomic cushion-top handle conforms to your hand and reduces impact shock
- Rubber surrounds on the handle allow resting the cane against walls or tables
- Height-adjustable design with 5 settings, suitable for users 152–198 cm tall
- Sturdy folding mechanism with secure strap for compact storage
- Flex-grip tip provides superior traction and smooth contact with walking surfaces
- Treaded base design improves stability and reduces slip risk
- Lightweight and strong construction for everyday mobility support

SPECIFICATIONS

- **Height Adjustment Range:** 32" – 37" (812 – 940 mm)
- **Height Settings:** 5
- **User Height Compatibility:** 152 – 198 cm
- **Weight Capacity:** 113 kg (250 lb)
- **Handle Type:** Ergonomic cushion-top
- **Tip Type:** Flex-grip with tread design
- **Colour:** Blue

- **Foldable:** Yes, with secure strap for storage