




## **WALKER WHEELCHAIR TRANSIT - AIRGO FUSION 2 IN 1**

**\$809.00**



**SKU: 700-935**

### **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Airgo Fusion Walker Wheelchair was designed to allow you to remain independent and live your life to the fullest. Its lightweight design and unparalleled stability provide the confidence to walk safely again. And if you ever find that you've ventured out too far, the Airgo Fusion converts easily to a transport chair so that a loved one can help bring you back home in comfort and style. Find out more about why so many people love their Airgo Fusion!

### **WALKER WHEELCHAIR FEATURES**

- Lightweight side-folding rollator and transport chair in one design
- Folds sideways easily with just a pull of the strap!
- Patent-pending Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest your arms in either mode
- Smooth bearings and an improved centre of gravity for easy pushing
- Unique ergonomic handgrips fit naturally for easier brake access
- Features our new senior-friendly tool-free brake adjustment knob
- The lightweight frame (17.1 lb net) designed for superior stability
- The reflective soft basket can be removed and is machine washable
- 8" (20 cm) Comfort-Grip wheels are ideal for indoor & outdoor use
- Net weight (without accessories): 17.1 lb (7.76 kg)

Buy this Airgo Fusion 2 in 1 Walker Transit Wheelchair online at MediMart using Afterpay or [NDIS method](#) today!

## **SPECIFICATIONS**

- Suitable for user height: 153 – 188cm
- Handle height adjustment: 76 – 93cm
- Seat height: 51cm
- Seat depth: 30cm
- Width between handles: 50cm
- Overall width: 69cm
- Distance between front & back wheels: 72cm
- Weight: 8.5kg
- Max. user weight: 136kg