


WALKER WHEELCHAIR TRANSIT - AIRGO FUSION 2 IN 1



SKU: 700-935



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Airgo Fusion Walker Wheelchair was designed to allow you to remain independent and live your life to the fullest. Its lightweight design and unparalleled stability provide the confidence to walk safely again. And if you ever find that you've ventured out too far, the Airgo Fusion converts easily to a transport chair so that a loved one can help bring you back home in comfort and style. Find out more about why so many people love their Airgo Fusion!

WALKER WHEELCHAIR FEATURES

- Lightweight side-folding rollator and transport chair in one design
- Folds sideways easily with just a pull of the strap!
- Patent-pending Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest your arms in either mode
- Smooth bearings and an improved centre of gravity for easy pushing
- Unique ergonomic handgrips fit naturally for easier brake access
- Features our new senior-friendly tool-free brake adjustment knob
- The lightweight frame (17.1 lb net) designed for superior stability
- The reflective soft basket can be removed and is machine washable
- 8" (20 cm) Comfort-Grip wheels are ideal for indoor & outdoor use
- Net weight (without accessories): 17.1 lb (7.76 kg)

Buy this Airgo Fusion 2 in 1 Walker Transit Wheelchair online at MediMart using Afterpay or [NDIS method](#) today!

SPECIFICATIONS

- Suitable for user height: 153 – 188cm
- Handle height adjustment: 76 – 93cm
- Seat height: 51cm
- Seat depth: 30cm
- Width between handles: 50cm
- Overall width: 69cm
- Distance between front & back wheels: 72cm
- Weight: 8.5kg
- Max. user weight: 136kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au