



AIRGO MINIQUAD ULTRA-STABLE CANE TIP

~~\$31.00~~ Original price was: \$31.00. \$27.90 Current price is: \$27.90.



SKU: 780-030

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Airgo Mini Quad Cane Tip** is a replacement base designed to provide enhanced stability and balance for your walking cane. Ideal for users who want the benefits of a quad base without the added bulk, this lightweight and compact tip helps keep your cane upright when not in use and ensures smoother, more secure movement on various surfaces.

FEATURES

- Compact MiniQuad design offers superior stability while remaining lightweight
- 54% lighter and 80% smaller than traditional large-base quad tips
- Self-centring tip maintains even contact with the ground for improved balance
- Keeps cane upright when not in use for convenience and safety
- Durable rubber construction for secure traction on most surfaces
- Compatible with most standard walking canes
- Quick and easy to attach or replace

SPECIFICATIONS

- **Tip Type:** MiniQuad – self-centring, stable quad design
- **Material:** Reinforced non-slip rubber
- **Compatibility:** Fits most 19 mm ($\frac{3}{4}$ ") cane shafts
- **Base Size:** Compact quad footprint
- **Colour:** Black
- **Use:** Replacement tip for canes requiring added stability