



ALIVIO DONATELLO PETITE PLUS

\$4,957.00



SKU: KA555-PL

PRODUCT DESCRIPTION

Purpose & Benefits of the Alivio Donatello Petite Plus

The Alivio Donatello Petite Plus is a state-of-the-art mobility and comfort aid, specifically designed for individuals with smaller frames who demand premium support and luxurious features. This advanced chair is perfect for both home and care environments, offering unparalleled customization and comfort to enhance your daily life.

Key Features

- **5 Independent Motors:** The chair is equipped with five powerful motors, allowing for infinite positioning options. You can easily adjust the headrest, backrest, lumbar support, footrest, and tilt position to achieve the perfect seating or reclining posture.
- **Zero Gravity Tilt-In-Space Function:** Experience the ultimate relaxation with the zero gravity feature, which creates a feeling of weightlessness, reducing pressure on the body and promoting a deep sense of calm.
- **Acupressure Massage:** The chair offers a soothing acupressure massage function with three adjustable pressure levels, designed to relieve stress and tension throughout your body.
- **Heating Function:** Integrated heating helps to reduce stress and pain while improving circulation, providing therapeutic warmth during use.
- **Lift and Tilt Function:** The unique lift and tilt feature ensures a smooth and effortless transition from sitting to standing, enhancing user independence and safety.
- **Premium Airline Head Support:** Designed for superior comfort, the head support offers extra cushioning for your head and neck, perfect for extended periods of use.
- **Durable and Stylish Fabric:** Upholstered in ultra-durable and stylish fabric, the Donatello Petite Plus not only looks good but is built to last.

- **Soft-Touch Hand Controller:** The intuitive hand controller features a backlit design for easy use in any lighting condition and includes an integrated USB charging point to keep your devices powered.
- **Retractable Drinks Holder:** A convenient drinks holder is built into the left side of the chair, retracting out of sight when not in use.
- **Wireless Phone Charger:** Located on the right side, this retractable wireless charger also includes controls for the heating and massage functions, making it easy to stay connected and relaxed.

Specifications

Overall Width	1010 mm
Overall Depth	1060 mm
Overall Height	1090 mm
Seat Width	520 mm
Seat Depth	520 mm
Floor To Seat Height	480 mm
Backrest Height	690 mm
Product Weight	79 kg
Max User Weight	158 kg
Brand	Alivio

The Alivio Donatello Petite Plus is the perfect choice for individuals with smaller body frames who seek a high-quality, comfortable, and reliable mobility aid. Its compact design and advanced features make it an excellent solution for promoting independence and comfort in various care environments.