



ALIVIO DONATELLO PLUS



SKU: KA554-PL

PRODUCT DESCRIPTION

Purpose & Benefits of the Alivio Donatello Plus

The Alivio Donatello Plus is a premium mobility and comfort solution, meticulously designed for those who demand advanced support and luxurious features in their daily lives. Whether in a home or care facility, this chair offers a perfect blend of functionality and comfort, making it an ideal choice for enhancing user independence and well-being.

Key Features

- **5 Independent Motors:** Experience unparalleled customization with five powerful motors that offer infinite positioning options. Easily adjust the headrest, backrest, lumbar support, footrest, and tilt to find your ideal seating or reclining posture.
- **Zero Gravity Tilt-In-Space Function:** Enjoy the sensation of weightlessness with the zero gravity feature, which minimizes pressure on your body and promotes deep relaxation.
- **Acupressure Massage:** The chair includes an acupressure massage function with three pressure levels, designed to relieve stress and tension, providing therapeutic comfort right at your fingertips.
- **Heating Function:** Integrated heating offers soothing warmth that helps reduce stress, alleviate pain, and improve circulation, making it perfect for those in need of therapeutic benefits.
- **Lift and Tilt Function:** The unique lift and tilt feature ensures a smooth and effortless transition between sitting and standing, enhancing both safety and independence.
- **Premium Airline Head Support:** Additional head and neck support is provided by the premium airline-style headrest, ensuring comfort during extended use.
- **Durable and Stylish Fabric:** The chair is upholstered in ultra-durable, stylish fabric that not

only enhances its appearance but also ensures longevity and ease of maintenance.

- **Soft-Touch Hand Controller:** The intuitive hand controller features a backlit design for easy use in low light conditions, and includes an integrated USB charging point, keeping your devices powered and within reach.
- **Retractable Drinks Holder:** Conveniently located on the left side, the drinks holder retracts when not in use, maintaining the chair's sleek design.
- **Wireless Phone Charger:** The right side of the chair features a retractable wireless charger, along with controls for the heating and massage functions, offering convenience and comfort at your fingertips.

Specifications

Total Height	1090 mm
Total Width	1010 mm
Total Depth	1060 mm
Seat Width	550 mm
Seat Depth	530 mm
Backrest Height	690 mm
Floor to Seat Height	495 mm
Colour	Canyon Steel
Weight	79 kg
SWL	158 kg

The Alivio Donatello Plus is an exceptional choice for individuals seeking a versatile, comfortable, and reliable mobility aid. Its advanced features and durable design make it ideal for enhancing quality of life and promoting independence in various care environments.

Please Note:

This model features hidden cup holders and heating controls neatly built into the armrests. To access them, simply pull out the armrests — a compartment will be revealed for your drink and the heating controls. When not in use, the armrests appear standard, keeping the design neat and discreet.

[KCare Alivio User Manual](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au