




ALIVIO DONATELLO RECLINER LIFT CHAIR

\$3,049.00 - \$3,096.00 Price range: \$3,049.00 through \$3,096.00



SKU: KA554-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Alivio Donatello Recliner Lift Chair offers the peak of position customisation. With four motors the Donatello gives independent adjustment of the headrest, lumbar, back and footrest areas. This chair is available with a choice of backrest styles - standard support backrest, high support backrest or lateral support backrest.

The hand control has soft grip buttons and allows for 6 pre-defined positions, allowing you to rise to a full standing position or a full recline, all with the touch of a button.

4 Independent motors for a wide range of position combinations and support:

Power Headrest - The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.

Power Lumbar - Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.

Power Backrest and Power Footrest - The dual lift system, allows you to individually adjust the backrest and the footrest for the position that best fits the clinical needs of the client.

Easy to use hand control

- Soft grip buttons

- 6 pre-defined positions
- Choose “full recline” or “full stand” with the touch of a button

[Click here for to view the full Alivio Donatello Recliner Brochure.](#)

The Alivio Donatello Recliner Lift Chair is a great choice for people who are looking for a comfortable, supportive, and easy-to-use lift chair. It is perfect for people with disabilities or mobility issues, as well as anyone who wants a chair that can provide maximum comfort and relaxation.

OPTIONS

The Donatello has three configurable backrest options to suit any need

FEATURES

- Dual lift function
- Lay-flat recliner
- Footrest extension
- All close & all lift buttons
- Hand controller
- Ultra-durable fabrics
- Lumbar support to accommodate sitting posture
- Standard support backrest, high support backrest or lateral (side) support backrest options

SPECIFICATIONS

| Specification | Standard Size | Petite Size |
|-------------------|---------------|-------------|
| Overall Depth | 1016mm | 978mm |
| Overall Width | 890mm | 851mm |
| Overall Height | 1092mm | 1041mm |
| Seat Depth | 520mm | 457mm |
| Seat Width | 520mm | 457mm |
| Seat Height | 508mm | 457mm |
| Reclined Length | 1830mm | 1830mm |
| Product Weight | 74kg | 70kg |
| Max User Weight | 158kg | 158kg |
| Safe Working Load | N/A | 158kg |