



ALIVIO MICHELANGELO RECLINER LIFT CHAIR



SKU: KA55-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF ALIVIO MICHELANGELO RECLINER LIFT CHAIR

The Alivio Michelangelo Recliner Lift Chair offers the peak of position customisation. With four motors, the Michelangelo gives independent adjustment of the headrest, lumbar, back, and footrest areas. The hand control has soft grip buttons and allows for 6 pre-defined positions, allowing you to rise to a full standing position or a full recline, all with the touch of a button.

The Alivio Michelangelo Recliner Lift Chair also features a power headrest and power lumbar support. The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back. Individually adjustable lumbar support can assist in maintaining proper sitting posture and relieving tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems. With the dual lift system, it allows you to individually adjust the backrest and the footrest for the position that best fits your clinical needs.

FEATURES

- **Power Lift and Recline Mechanism:** Smooth and quiet operation allows for effortless adjustment to any position, from fully reclined to a gentle lift that aids in standing up.
- **Luxurious Upholstery:** Wrapped in high-quality, soft fabric that is both durable and easy to clean, ensuring long-lasting comfort and appeal.
- **Ergonomic Design:** Thoughtfully designed with user comfort in mind, featuring plush padding, a supportive backrest, and a spacious seating area.
- **Built-in Remote Control:** Conveniently located and easy-to-use remote for adjusting the chair's positions and activating the lift function without any strain.

- **Safety Features:** Equipped with emergency battery backup and anti-tipping design, ensuring users can safely enjoy the chair's full range of features even during a power outage.
- **Multiple Reclining Positions:** Allows users to find their perfect comfort spot for reading, watching TV, or taking a restful nap.

EASY TO USE HAND CONTROL

- Soft grip buttons
- 6 pre-defined positions
- Choose "full recline" or "full stand" with the touch of a button.

SPECIFICATIONS

- **Material:** Premium fabric upholstery with metal and wood frame
- **Power Source:** Standard electrical outlet, includes emergency battery backup
- **Size Available:** Standard, Petite

Dimension	Standard	Petite
Overall Length	1080 mm	-
Max Overall Length	1830 mm	1830 mm
Overall Width	890 mm	812 mm
Overall Depth	1080 mm	978 mm
Overall Height	1067 mm	1016 mm
Seat Width	559 mm	457 mm
Seat Depth	546 mm	457 mm
Max Seat Height	470 mm	457 mm
Product Weight	70 kg	-
Max User Weight	158 kg	136 kg

[Click here for specifications and further information.](#)

[KCare Alivio User Manual](#)

VIDEO

https://youtu.be/_msquBJ7r6s

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au