



ALLCARE RESISTANT BANDS - SIZE & COLOUR OPTIONS



SKU: ACBAND5

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Allband Exercise Band is used in the rehabilitation and strengthening of joints in the musculature system. Allband Exercise Band can be used with ROM exercises to provide resistance and a resultant increase in strength. Use of Allband Exercise Band can be throughout the rehabilitation process as it comes in varying grades of resistance, ranging from easy to super heavy. **The resistance is marked by colour of the band, light colours = light resistance, and dark colours = heavy resistance.**

FEATURES

- 5.5 metre or 1.5 metre length
- Steady 30% increase from colour to colour
- To use, simply cut a length from the roll and tie a knot to form a loop. The knot can be undone easily to change the length or size of the loop. Allband Exercise Band may also be used as a single length
- Available in 6 colour coded resistances, providing a wide range of treatments
- Made under the international quality standard ISO 9002
- Supplied in colour coded boxes, with perforated dispensing slot for easy access and storage

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au