



## ALLCARE RESISTANT BANDS - SIZE & COLOUR OPTIONS

\$16.00



**SKU:** ACBAND5

## **PRODUCT DESCRIPTION**

## **PURPOSE AND BENEFITS**

Allband Exercise Band is used in the rehabilitation and strengthening of joints in the musculature system. Allband Exercise Band can be used with ROM exercises to provide resistance and a resultant increase in strength. Use of Allband Exercise Band can be throughout the rehabilitation process as it comes in varying grades of resistance, ranging from easy to super heavy. **The resistance is marked by colour of the band, light colours = light resistance, and dark colours = heavy resistance.** 

## **FEATURES**

- 5.5 metre or 1.5 metre length
- Steady 30% increase from colour to colour
- To use, simply cut a length from the roll and tie a knot to form a loop. The knot can be undone easily to change the length or size of the loop. Allband Exercise Band may also be used as a single length
- Available in 6 colour coded resistances, providing a wide range of treatments
- Made under the international quality standard ISO 9002
- Supplied in colour coded boxes, with perforated dispensing slot for easy access and storage