



# ALLEGRO ALTO 200 MOBILE HOIST (EXCL. SPREADER BAR)

\$4,416.00 Original price was: \$4,416.00.<u>\$3,974.00</u>Current price is: \$3,974.00.



SKU: ALTO01ELS

## PRODUCT DESCRIPTION

## **PURPOSE & BENEFITS**

The Allegro ALTO 200 Mobile Hoist is an innovative patient-lifting solution that prioritizes safety, functionality, and convenience. Engineered to make patient transfers effortless and secure, it's an invaluable asset for caregivers in both home and clinical settings.

#### **FEATURES**

- Lifting Capacity: Designed to safely lift and transport patients up to 200kg.
- **Easy Maneuverability:** Fitted with premium quality, lockable castor wheels for smooth transitions and steering.
- Ergonomic Handset: Intuitive controls for caregivers, ensuring precision during lifts.
- Emergency Stop & Lower: Built-in safety mechanisms allow for an immediate stop or controlled lowering in emergencies.
- Adjustable Boom Height: Offers flexible lifting range, catering to various transfer scenarios.
- Compact Design: Ideal for use in confined spaces or rooms with limited maneuverability.
- **Rechargeable Battery System:** Ensures uninterrupted operations, with visual and auditory low battery indicators.

#### **BENEFITS**

- **Safety First:** Equipped with various safety features, like anti-pinch and anti-swing, ensuring utmost patient security.
- User-Friendly: Intuitive design and simple controls make it easy for caregivers to operate, reducing training time.
- Versatile Use: Suitable for various transfer situations, including bed-to-chair, floor lifts, and



bathroom transfers.

• **Durable Construction:** Made with high-quality materials that guarantee longevity and sturdiness.

# **APPLICATIONS**

The Allegro ALTO 200 Mobile Hoist is suitable for:

- Hospitals and clinical environments where patient mobility is paramount.
- Rehabilitation centres and physiotherapy clinics.
- Assisted living and elderly care facilities.
- Home use for patients with severe mobility restrictions.

#### **SPECIFICATIONS**

Safe Working Load: 200kg

• Base Length: 1150mm

• External Base Width Closed: 700mm

• External Base Width Open: 1130mm

• Internal Base Width Closed: 600mm

• Internal Base Width Open: 1030mm

• Base Width Adjustment: Electric

• Under Bed Clearance: 130mm

• Minimum Height Yoke: 300mm

• Maximum Height Yoke: 1920mm

• Average Lift Range\*: 1620mm

• Lifter Weight: 46kg

• Minimum Storage Height: 1350mm

• Under Base Height: 55mm

• Maximum Hoist Reach: 1060mm

• Minimum Hoist Reach: 610mm

• Hoist Reach at Maximum Height: 515mm

Mass of Heaviest Part: 24kg (Mast assembly)

### see options below for optional extra:

- Allegro Standard Yoke with Integrated Weigh Scale Spreader Bar
- Allegro Tri Hook Yoke Spreader Bar
- Allegro 4 Point Yoke Spreader Bar
- Allegro Standard Yoke Spreader Bar