




ALLEGRO ALTO 200 MOBILE HOIST (EXCL. SPREADER BAR)



SKU: ALTO01ELS

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Allegro ALTO 200 Mobile Hoist is an innovative patient-lifting solution that prioritizes safety, functionality, and convenience. Engineered to make patient transfers effortless and secure, it's an invaluable asset for caregivers in both home and clinical settings.

FEATURES

- **Lifting Capacity:** Designed to safely lift and transport patients up to 200kg.
- **Easy Maneuverability:** Fitted with premium quality, lockable castor wheels for smooth transitions and steering.
- **Ergonomic Handset:** Intuitive controls for caregivers, ensuring precision during lifts.
- **Emergency Stop & Lower:** Built-in safety mechanisms allow for an immediate stop or controlled lowering in emergencies.
- **Adjustable Boom Height:** Offers flexible lifting range, catering to various transfer scenarios.
- **Compact Design:** Ideal for use in confined spaces or rooms with limited maneuverability.
- **Rechargeable Battery System:** Ensures uninterrupted operations, with visual and auditory low battery indicators.

BENEFITS

- **Safety First:** Equipped with various safety features, like anti-pinch and anti-swing, ensuring utmost patient security.
- **User-Friendly:** Intuitive design and simple controls make it easy for caregivers to operate, reducing training time.
- **Versatile Use:** Suitable for various transfer situations, including bed-to-chair, floor lifts, and

bathroom transfers.

- **Durable Construction:** Made with high-quality materials that guarantee longevity and sturdiness.

APPLICATIONS

The Allegro ALTO 200 Mobile Hoist is suitable for:

- Hospitals and clinical environments where patient mobility is paramount.
- Rehabilitation centres and physiotherapy clinics.
- Assisted living and elderly care facilities.
- Home use for patients with severe mobility restrictions.

SPECIFICATIONS

- **Safe Working Load:** 200kg
- **Base Length:** 1150mm
- **External Base Width Closed:** 700mm
- **External Base Width Open:** 1130mm
- **Internal Base Width Closed:** 600mm
- **Internal Base Width Open:** 1030mm
- **Base Width Adjustment:** Electric
- **Under Bed Clearance:** 130mm
- **Minimum Height Yoke:** 300mm
- **Maximum Height Yoke:** 1920mm
- **Average Lift Range*:** 1620mm
- **Lifter Weight:** 46kg
- **Minimum Storage Height:** 1350mm
- **Under Base Height:** 55mm
- **Maximum Hoist Reach:** 1060mm
- **Minimum Hoist Reach:** 610mm
- **Hoist Reach at Maximum Height:** 515mm
- **Mass of Heaviest Part:** 24kg (Mast assembly)

see options below for optional extra:

- [Allegro Standard Yoke with Integrated Weigh Scale Spreader Bar](#)
- [Allegro Tri Hook Yoke Spreader Bar](#)
- [Allegro 4 Point Yoke Spreader Bar](#)
- [Allegro Standard Yoke Spreader Bar](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au