



## ALLEGRO AMPUTEE SLING

**\$345.00 - \$469.00** Price range: \$345.00 through \$469.00



**SKU:** A0-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Allegro Amputee Slings are recommended for lifting patients with high-level amputations. These slings incorporate a head support and a commode aperture, providing a safe and secure lift while maintaining full functionality. Designed to be simple and safe to use for both home and institutional patient lifting, these slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas. These slings utilise webbing attachment loops, which are colour coded for sizing to allow easy sling selection and identification. Webbing loops of varying lengths are provided at each connection point to allow a layback posture for the patient or, by using the shortest loops, a higher lift.

### FEATURES

- Sling incorporates a head support and commode aperture
- Recommended for lifting patients with high-level amputations
- Safe Working Load: 300kg

## SIZING