



ALLEGRO AMPUTEE SLING



SKU: A0-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Allegro Amputee Slings are recommended for lifting patients with high-level amputations. These slings incorporate a head support and a commode aperture, providing a safe and secure lift while maintaining full functionality. Designed to be simple and safe to use for both home and institutional patient lifting, these slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas. These slings utilise webbing attachment loops, which are colour coded for sizing to allow easy sling selection and identification. Webbing loops of varying lengths are provided at each connection point to allow a layback posture for the patient or, by using the shortest loops, a higher lift.

FEATURES

- Sling incorporates a head support and commode aperture
- Recommended for lifting patients with high-level amputations
- Safe Working Load: 300kg

SIZING

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au