



## **ALLEGRO DELUXE LEG SLING**

**\$401.00 - \$487.00** Price range: \$401.00 through \$487.00



**SKU:** DL0-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Deluxe Leg Slings are designed to provide maximum support to hips & thighs. These slings are simple & safe to use for both home & institutional patient lifting.

The slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas. These slings utilise webbing attachment loops, which are colour coded for sizing to allow easy sling selection & identification. Webbing loops of varying lengths are provided at each connection point to allow a lay back posture for the patient or, by using the shortest loops, a higher lift.

### **FEATURES**

- Provides maximum support to hips and thighs.
- Recommended for patients with compromised hip joints and amputees (subject to assessment).
- Safe Working Load: 300kg

### **SIZING**