



ALLEGRO DELUXE LEG SLING



SKU: DL0-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Deluxe Leg Slings are designed to provide maximum support to hips & thighs. These slings are simple & safe to use for both home & institutional patient lifting.

The slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas. These slings utilise webbing attachment loops, which are colour coded for sizing to allow easy sling selection & identification. Webbing loops of varying lengths are provided at each connection point to allow a lay back posture for the patient or, by using the shortest loops, a higher lift.

FEATURES

- Provides maximum support to hips and thighs.
- Recommended for patients with compromised hip joints and amputees (subject to assessment).
- Safe Working Load: 300kg

SIZING

For further information go to medimart.com.au or give us a call: 1300 633 063

SA



Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au