



## ALLEGRO FORTE 320 BARIATRIC LIFTER

\$5,717.00



SKU: FOR01LW

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Forte 320 Bariatric Hoist is a robust and reliable patient lifting device designed to safely and comfortably handle the specific needs of bariatric patients. This hoist is an essential tool in healthcare settings, ensuring the dignified and secure transfer of individuals with higher body weight.

**Bariatric Yoke Spreader Bar 4 Point sold separately. [Click Here](#)**

### FEATURES

- **High Weight Capacity:** Engineered to support substantial weight, the Forte 320 is ideal for bariatric patients, providing a safe lifting capacity that meets the demands of heavier weight ranges.
- **Four-Point Yoke System:** This feature offers enhanced stability and even weight distribution, reducing discomfort for the patient during lifting and transfers.
- **Sturdy and Durable Design:** Constructed with high-quality materials, the hoist is built for durability and long-term use, ensuring reliability in a professional healthcare environment.
- **Ease of Manoeuvrability:** Despite its robust construction, the hoist is designed for ease of movement, with smooth-rolling casters that make it easy to navigate across various floor surfaces.
- **Simple and Intuitive Controls:** The hoist is equipped with user-friendly controls, making it easy for caregivers to operate, and ensuring safe and precise handling.
- **Emergency Safety Features:** Includes emergency lowering and stop functions for additional safety, providing peace of mind in case of power loss or other unexpected situations.
- **Comfortable Patient Sling Options:** Compatible with a range of sling types and sizes, allowing for a customised fit and ensuring the comfort and security of the patient during

transfers.

## **SPECIFICATIONS**

- **Safe Working Load:** 320 kg
- **Base Length:** 1420 mm
- **External Base Width Closed:** 715 mm
- **External Base Width Open:** 1200 mm
- **Internal Base Width Closed:** 600 mm
- **Internal Base Width Open:** 1080 mm
- **Base Width Adjustment:** Electric
- **Low Profile Leg:** Twin 75 mm Castors
- **Under Bed Clearance:** 111 mm
- **Minimum Yoke Height:** 245 mm
- **Maximum Yoke Height:** 1940 mm
- **Overall Height (Boom Raised):** 2243 mm
- **Lifter Weight:** 63.0 kg
- **Minimum Storage Height:** 1445 mm
- **Under Base Height:** 15 mm
- **Maximum Hoist Reach:** 810 mm
- **Mass of Heaviest Part (Chassis):** 37 kg

The Forte 320 Bariatric Hoist is a vital piece of equipment in the field of bariatric care, offering a combination of strength, safety, and ease of use. It is designed to enhance the quality of care for patients with higher weight needs, ensuring their comfort and dignity during transfers.