



ALLEGRO GENERAL PURPOSE SLING

\$263.00 - \$361.00 Price range: \$263.00 through \$361.00



SKU: GP0-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Allegro General Purpose Slings are designed to be simple and safe to use for both home and institutional patient lifting. The slings are manufactured from polyester fabric or mesh (suitable for bathing) with internal padding in pressure areas and utilise webbing attachment loops. The webbing loops are colour coded for sizing to allow easy sling selection and identification. Webbing loops of varying lengths are provided at each connection point to allow a lay back posture for the patient or, by using the shortest loops, a higher lift.

FEATURES

Six sizes in fabric or mesh.

- Provides easy fitting and good
- Recommended for general lifting
- Safe Working Load: 300kg

SIZING